

# Sterilization Container Systems

---

Secure Today  
Elevate Tomorrow

---


 Different  
Colors





Victoris retains the right to modify and enhance the products showcased in this catalogue as part of ongoing development efforts. All provided weight and dimensional information are approximate and may vary from the actual products.

Copyright © Victoris Co., Inc. All Rights Reserved  
VH-CC-K001 Rev.00



## Table Of Contents

---

|           |                          |           |
|-----------|--------------------------|-----------|
| <b>06</b> | INTRODUCTION             | <b>15</b> |
| <b>16</b> | STANDARD MODELS          | <b>35</b> |
| <b>36</b> | REUSABLE MODELS          | <b>57</b> |
| <b>58</b> | STERILIZATION BASKETS    | <b>67</b> |
| <b>68</b> | ACCESSORIES & COMPONENTS | <b>69</b> |
| <b>70</b> | TROLLEYS & TABLES        | <b>71</b> |
| <b>72</b> | INDEX                    | <b>73</b> |

## The Full Story

Victoris is more than a healthcare supply company; it's a testament to resilience, integrity, and a commitment to human values. Founded by two industry professionals who shared a vision for a workplace that prioritizes people over profit, our story began in the halls of an international medical device company. In the face of adversity, to create a company that not only thrives in the industry but serves as a beacon of compassion and ethical leadership. Victoris was born from the belief that success is measured not just in financial gains but in the positive impact a company can have on its employees and the community it serves.



### Our Vision

At Victoris, our vision is to be a leader in healthcare solutions, recognized not only for the quality of our products but for the unwavering commitment to the well-being of our employees and the communities we serve. We envision a future where healthcare is not only about cutting-edge technology but also about compassionate and ethical practices, setting a benchmark for the industry.

### Our Mission

Our mission at Victoris is twofold. Firstly, to provide exceptional healthcare products that enhance the lives of individuals undergoing medical procedures, for humans or pets. Secondly, to pioneer a corporate culture that prioritizes the welfare of our team members, recognizing them as the driving force behind our success. We are dedicated to fostering an environment that values diversity, upholds human rights, and serves as a beacon for responsible and ethical business practices.

## Precision in Form, Compliance in Function

### R&D

Experience the innovation driving Victoris with our dedicated Research and Development (R&D) division. Our R&D team is committed to pushing the boundaries of healthcare technology, continually developing breakthrough solutions to meet the evolving needs of the industry.

Utilizing cutting-edge technologies and advanced methodologies, our focus on rigorous testing and constant improvement ensures that each product meets the highest standards of quality and reliability. Trust in our relentless pursuit of excellence to elevate your operations and achieve unparalleled performance.

### Production

Discover the excellence of Victoris production processes, designed to ensure precision, efficiency, and unparalleled quality. Our advanced production techniques are meticulously developed to deliver consistent results, ensuring that every product meets the highest standards of safety and reliability.

Spanning a closed production area of 2000 square meters, we leverage innovative technologies and streamlined workflows to minimize downtime and enhance productivity. Trust in our expertise to maintain the utmost quality and reliability in every aspect of your operations.

### Regulatory / Quality

At Victoris, our unwavering commitment to quality, safety, and compliance, including adherence to regulatory bodies such as CE, FDA, and ISO, is the driving force behind everything we do. Our Regulatory Affairs and Quality Assurance department, at the heart of our operations, ensures that our products consistently meet the highest industry standards and comply with all relevant guidelines.

We understand that in the medical industry, quality is non-negotiable, which is why our RA/QA team is dedicated to upholding the most rigorous standards of product quality, safety, and effectiveness, leaving no stone unturned to meet or exceed the stringent requirements of the healthcare and veterinary sectors.

Navigating the complex landscape of regulatory compliance is essential, and our experienced regulatory affairs professionals stay up to date with the latest regulations, working diligently to secure the necessary approvals, certifications, and registrations from CE, FDA, ISO, and other relevant authorities.

Furthermore, our commitment to Occupational Health and Safety Policy is integral to our operations, ensuring the well-being of our workforce and the quality of our products.

### Technical Competencies

- Machining
- Subtractive Manufacturing
- Deep Drawing
- Die & Mold Making
- Chemical Coating
- Hard Anodization
- Silicone Forming
- Laser Welding
- OEM - ODM - PLM



Special Stackability Design

2 mm Aluminium Bottom & Lid


OEM-PLM Marking Option



Place and Go Gasket

Shape / No Glue

Double-Button Extra-Secure Lock



Filter Retainer With Sealing & Spring System

Spring-Loaded Aluminium Lock for Better Grip and Sealing

Face Plate Made of Single Block of Aluminium

Extra-Strong Silicone Coated Ergonomic Handles



## Different Lids Colors

Our sterilization containers at Victoris come in eight different colors (silver, yellow, red, blue, green, black, bronze and purple) ensuring to easily distinguish between different sets or types of instruments at a glance.



### Advantages of Multiple Colors


- **Enhanced Organization:** Different colors can be used to categorize and organize instruments based on their usage, department, or level of sterilization, reducing the risk of cross-contamination and improving efficiency.
- **Easy Identification:** Color-coding helps in quickly identifying the contents of each container, saving time and reducing errors in busy environments.
- **Improved Workflow:** With a clear visual system in place, staff can more easily manage inventory and ensure that the right instruments are always in the right place.

## Customized Solutions & Tailored to Your Needs

Victoris sterilization containers are manufactured in 8 main sizes to accommodate various surgical needs: standard sizes (1/1 Full Size, 3/4 Quarter Size, 1/2 Half Size, Mini, Dental and 1/2 Dental Size) are ideal for common instrument sets, while the Wide Body size suits specialized instruments like Bookwalter retractors, and the Long Body size is designed for extended instruments such as robotic arms.

The sterilization containers are crafted in five primary heights to suit a variety of surgical instruments. Additionally, custom sizes or heights can be tailored to meet specific requirements.

Standard 5 Different Heights For Standard & Reusable Models



Standard 8 Different Size Types For Standard & Reusable Models

| Type  | Packaging according to EN 868-8   |                    | Steam sterilization procedure according to EN 285 |            |                   |                 |                 |
|---|---|--------------------|---|------------|-------------------|-----------------|-----------------|
|   | Description   | Symbolic Packaging | Gravity   | Pre-Vacuum | Fractional Vacuum | Fractional Flow | Steam Injection |
| Sterilization Container with filters  | Textiles in container with filter in lid and box                                |                    | +   | +          | +                 | +               | +               |
|   |   |                    | O   | O          | +                 | +               | +               |
|   | Textiles in container with filter in lid, non-perforated box                    |                    | -   | -          | +                 | +               | +               |
|   |   |                    | -   | -          | +                 | +               | O               |
|   | Textiles in container with filter in box, non-perforated lid                    |                    | O   | +          | +                 | +               | +               |
|   |   |                    | -   | O          | +                 | +               | +               |
|   | Instruments in wire basket, in sterilization filter in lid and box              |                    | +3  | +          | +                 | +               | +               |
| Instruments in wire basket, in sterilization container filter in lid and non-perforated box     |   | -                  | +2  | +2         | +2                | +2              |                 |
| Instruments in wire basket, in sterilization container filter in box and non-perforated lid     |   | O                  | +   | +          | +                 | +               |                 |
| Sterilization Container with valves   | Textiles in sterilization container with valves in lid and box                  |                    | -   | O          | +                 | +               | O               |
|   |   |                    | -   | -          | +                 | +               | O               |
|   | Textiles in sterilization container with valves in lid and no valves on the box |                    | -   | -          | +                 | +               | -               |
|   |   |                    | -   | -          | +                 | +               | -               |
|   | Instruments in wire basket in sterilization container with valve in lid and box |                    | -   | +2         | +2                | +2              | +2              |
| Instruments in wire basket, in sterilization container with valve in lid and non-perforated box |   | -                  | +2  | +2         | +2                | +2              |                 |

\*To prevent errors during service, operation, or loading of the autoclaves, consult above table.  
 \*Use only sterilization containers that guarantee proper evaporation or discharge of condensate.  
 \*Keep in mind that gravitation-based procedures, which lack a dedicated drying process, do not ensure sufficient after-drying effects.

# STANDARD & REUSABLE MODELS

SECURE TODAY

ELEVATE TOMORROW



## WHY STERILIZATION CONTAINERS?



- Consistent Processing
- Reduced Risk of Contamination
- Physical Instrument Protection
- Optimal Set Organization
- Reusability / Sustainability
- Reduced Waste
- Streamlined Process
- Easy Identification and Access
- Adherence to Standards
- Minimized Handling

# Today's Safety. Tomorrow's Efficiency



## Elevated Sterilization Solutions: For Your Operational Excellence

Discover our 2mm special alloy aluminum container systems that provide superior protection for your surgical instruments. Available in 8 different color palettes, they facilitate organization and quick identification. This solution supports your operational excellence by combining reliability, an eco-friendly approach, and long-term cost advantages.



## Standard Model

1/1 Full Size

 Different Lids Colors



With Safety Lid



With Perforated Lid



Non Perforated



**Code**      **Outside Dimension**  
**Inside Dimension**  
**L x W x H**

K1-10-110    595 x 270 x 95 mm  
544 x 257 x 92 mm

K1-10-135    595 x 270 x 115 mm  
544 x 257 x 113 mm

K1-10-150    595 x 270 x 140 mm  
544 x 257 x 137 mm

K1-10-200    595 x 270 x 190 mm  
544 x 257 x 187 mm

K1-10-260    595 x 270 x 250 mm  
544 x 257 x 247 mm



Perforated



**Code**      **Outside Dimension**  
**Inside Dimension**  
**L x W x H**

KP1-10-110    595 x 270 x 95 mm  
544 x 257 x 80 mm

KP1-10-135    595 x 270 x 115 mm  
544 x 257 x 101 mm

KP1-10-150    595 x 270 x 140 mm  
544 x 257 x 125 mm

KP1-10-200    595 x 270 x 190 mm  
544 x 257 x 175 mm

KP1-10-260    595 x 270 x 250 mm  
544 x 257 x 235 mm



Lids



**Colors**    **Code**

Grey    K1-10-1  
Yellow    K1-10-2  
Red    K1-10-3  
Blue    K1-10-4  
Green    K1-10-5  
Black    K1-10-6  
Purple    K1-10-7  
Bronze    K1-10-8

**Dimension**  
585 x 285 mm



Safety Lids




**Colors**    **Code**


Grey    KS1-10-1  
Yellow    KS1-10-2  
Red    KS1-10-3  
Blue    KS1-10-4  
Green    KS1-10-5  
Black    KS1-10-6  
Purple    KS1-10-7  
Bronze    KS1-10-8

**Dimension**  
605 x 287 mm

### Wire Basket



Mesh: 4x4 mm  
Wire: ø 1 mm



Box

| Code        | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B1-5425-30  | 540 x 255 x 30 mm              |
| B1-5425-50  | 540 x 255 x 50 mm              |
| B1-5425-70  | 540 x 255 x 70 mm              |
| B1-5425-100 | 540 x 255 x 100 mm             |




Lid

| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| B1-5425 | 510 x 240 x 5 mm               |

### Perforated Basket



Mesh: 4x4 mm  
Wire: ø 1 mm



Box (with 10 mm feet)

| Code        | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B2-5425-30  | 540 x 253 x 40 mm              |
| B2-5425-50  | 540 x 253 x 60 mm              |
| B2-5425-70  | 540 x 253 x 80 mm              |
| B2-5425-100 | 540 x 253 x 110 mm             |




Lid

| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| B2-5425 | 540 x 256 x 5 mm               |

### Hybrid Basket



Mesh: 4x4 mm  
Sides: ø 3 mm  
Wire: ø 1 mm



Box

| Code        | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B3-5425-50  | 540 x 255 x 50 mm              |
| B3-5425-70  | 540 x 255 x 70 mm              |
| B3-5425-100 | 540 x 255 x 100 mm             |



Lid

| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| B3-5425 | 540 x 255 x 5 mm               |

### Silicone Mat



| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| A7-5425 | 500 x 235 mm                   |


## Standard Model

1/1 Full Size Baskets

### Wire Basket



### Perforated Basket



### Hybrid Basket



## Standard Model

3/4 Quarter Size

 Different Lids Colors



With Safety Lid



With Perforated Lid



Non Perforated



**Code**      **Outside Dimension**  
**Inside Dimension**  
**L x W x H**

K1-20-110    470 x 270 x 95 mm  
430 x 260 x 92 mm

K1-20-135    470 x 270 x 115 mm  
430 x 260 x 101 mm

K1-20-150    470 x 270 x 140 mm  
430 x 260 x 137 mm

K1-20-200    470 x 270 x 184 mm  
430 x 260 x 178 mm



Perforated



**Code**      **Outside Dimension**  
**Inside Dimension**  
**L x W x H**

KP1-20-110    470 x 270 x 95 mm  
430 x 260 x 80 mm

KP1-20-135    470 x 270 x 115 mm  
430 x 260 x 101 mm

KP1-20-150    470 x 270 x 140 mm  
430 x 260 x 125 mm

KP1-20-200    470 x 270 x 184 mm  
430 x 260 x 166 mm



Lids



**Colors**    **Code**

Grey    K1-20-1  
Yellow    K1-20-2  
Red    K1-20-3  
Blue    K1-20-4  
Green    K1-20-5  
Black    K1-20-6  
Purple    K1-20-7  
Bronze    K1-20-8

**Dimension**  
465 x 285 mm



Safety Lids




**Colors**    **Code**


Grey    KS1-20-1  
Yellow    KS1-20-2  
Red    KS1-20-3  
Blue    KS1-20-4  
Green    KS1-20-5  
Black    KS1-20-6  
Purple    KS1-20-7  
Bronze    KS1-20-8

**Dimension**  
482 x 287 mm

### Wire Basket



Mesh: 4x4 mm  
Wire: ø 1 mm



Box

| Code        | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B1-4025-30  | 405 x 255 x 30 mm              |
| B1-4025-50  | 405 x 255 x 50 mm              |
| B1-4025-70  | 405 x 255 x 70 mm              |
| B1-4025-100 | 405 x 255 x 100 mm             |



Lid

| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| B1-4025 | 375 x 240 x 5 mm               |

### Perforated Basket



Mesh: 4x4 mm  
Wire: ø 1 mm



Box (with 10 mm feet)

| Code        | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B2-4025-30  | 405 x 252 x 40 mm              |
| B2-4025-50  | 405 x 252 x 60 mm              |
| B2-4025-70  | 405 x 252 x 80 mm              |
| B2-4025-100 | 405 x 252 x 110 mm             |




Lid

| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| B2-4025 | 409 x 256 x 5 mm               |

### Hybrid Basket



Mesh: 4x4 mm  
Sides: ø 3 mm  
Wire: ø 1 mm



Box

| Code        | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B3-4025-50  | 405 x 245 x 50 mm              |
| B3-4025-70  | 405 x 245 x 70 mm              |
| B3-4025-100 | 405 x 245 x 100 mm             |



Lid

| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| B3-4025 | 405 x 245 x 5 mm               |

### Silicone Mat



**Code**  
A7-4025  
**Outside Dimension**  
L x W x H  
370 x 235 mm

## Standard Model

3/4 Quarter Size Baskets

### Wire Basket



### Perforated Basket



### Hybrid Basket



## Standard Model

1/2 Half Size

 Different Lids Colors



With Safety Lid



With Perforated Lid



Non Perforated



**Code**      **Outside Dimension**  
**Inside Dimension**  
**L x W x H**

K1-30-110    300 x 270 x 95 mm  
255 x 255 x 92 mm

K1-30-135    300 x 270 x 115 mm  
255 x 255 x 113 mm

K1-30-150    300 x 270 x 140 mm  
255 x 255 x 137 mm

K1-30-200    300 x 270 x 190 mm  
255 x 255 x 187 mm

K1-30-260    300 x 270 x 250 mm  
255 x 255 x 247 mm



Perforated



**Code**      **Outside Dimension**  
**Inside Dimension**  
**L x W x H**

KP1-30-110    300 x 270 x 95 mm  
255 x 255 x 80 mm

KP1-30-135    300 x 270 x 115 mm  
255 x 255 x 101 mm

KP1-30-150    300 x 270 x 140 mm  
255 x 255 x 125 mm

KP1-30-200    300 x 270 x 190 mm  
255 x 255 x 175 mm

KP1-30-260    300 x 270 x 250 mm  
255 x 255 x 235 mm



Lids



**Colors**    **Code**

Grey    K1-30-1  
Yellow    K1-30-2  
Red    K1-30-3  
Blue    K1-30-4  
Green    K1-30-5  
Black    K1-30-6  
Purple    K1-30-7  
Bronze    K1-30-8

**Dimension**  
295 x 285 mm



Safety Lids




**Colors**    **Code**


Grey    KS1-30-1  
Yellow    KS1-30-2  
Red    KS1-30-3  
Blue    KS1-30-4  
Green    KS1-30-5  
Black    KS1-30-6  
Purple    KS1-30-7  
Bronze    KS1-30-8

**Dimension**  
315 x 287 mm

### Wire Basket



Mesh: 4x4 mm  
Wire: ø 1 mm



Box


| Code        | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B1-2525-30  | 255 x 245 x 30 mm              |
| B1-2525-50  | 255 x 245 x 50 mm              |
| B1-2525-70  | 255 x 245 x 70 mm              |
| B1-2525-100 | 255 x 245 x 100 mm             |




Lid

| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| B1-2525 | 240 x 215 x 5 mm               |

### Perforated Basket



Mesh: 4x4 mm  
Wire: ø 1 mm



Box (with 10 mm feet)

| Code        | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B2-2525-30  | 242 x 252 x 40 mm              |
| B2-2525-50  | 242 x 252 x 60 mm              |
| B2-2525-70  | 242 x 252 x 80 mm              |
| B2-2525-100 | 242 x 252 x 110 mm             |




Lid

| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| B2-2525 | 246 x 256 x 5 mm               |

### Hybrid Basket



Mesh: 4x4 mm  
Sides: ø 3 mm  
Wire: ø 1 mm



Box

| Code        | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B3-2525-50  | 255 x 245 x 50 mm              |
| B3-2525-70  | 255 x 245 x 70 mm              |
| B3-2525-100 | 255 x 245 x 100 mm             |



Lid

| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| B3-2525 | 245 x 245 x 5 mm               |

### Silicone Mat




| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| A7-2525 | 225 x 235 mm                   |


## Standard Model

1/2 Half Size Baskets

### Wire Basket



### Perforated Basket



### Hybrid Basket



## Standard Model

Long Size

 Different Lids Colors



With Perforated Lid



With Color Options



**Code**      **Outside Dimension**  
**Inside Dimension**  
**L x W x H**

K1-40-150    770 x 270 x 135 mm  
720 x 260 x 130 mm

K1-40-200    770 x 270 x 190 mm  
720 x 260 x 185 mm



**Code**      **Outside Dimension**  
**Inside Dimension**  
**L x W x H**

KP1-40-150    770 x 270 x 135 mm  
720 x 260 x 118 mm

KP1-40-200    770 x 270 x 190 mm  
720 x 260 x 173 mm



**Colors**    **Code**

Grey      K1-40-1

Yellow    K1-40-2

Red        K1-40-3

Blue       K1-40-4

Green     K1-40-5

Black      K1-40-6


Purple    K1-40-7

Bronze    K1-40-8

**Dimension**  
755 x 285 mm



### Wire Basket



**Code**      **Outside Dimension**  
**L x W x H**

B1-6625-50    665 x 250 x 50 mm  
B1-6625-70    665 x 250 x 70 mm  
B1-6625-100   665 x 250 x 100 mm

Mesh: 4x4 mm  
Wire: ø 1 mm



**Code**  
B1-6625

**Outside Dimension**  
**L x W x H**  
635 x 235 x 5 mm

### Perforated Basket



**Code**      **Outside Dimension**  
**L x W x H**

B2-6625-30    665 x 253 x 40 mm  
B2-6625-50    665 x 253 x 60 mm  
B2-6625-70    665 x 253 x 80 mm  
B2-6625-100   665 x 253 x 110 mm


Mesh: 4x4 mm  
Wire: ø 1 mm



**Code**  
B2-6625

**Outside Dimension**  
**L x W x H**  
668 x 256 x 5 mm

### Hybrid Basket



**Code**      **Outside Dimension**  
**L x W x H**

B3-6625-50    665 x 245 x 50 mm  
B3-6625-70    665 x 245 x 70 mm  
B3-6625-100   665 x 245 x 100 mm

Mesh: 4x4 mm  
Sides: ø 3 mm  
Wire: ø 1 mm



**Code**  
B3-6625

**Outside Dimension**  
**L x W x H**  
668 x 256 x 5 mm


## Standard Model

Long Size Baskets

### Wire Basket



### Perforated Basket



### Hybrid Basket



## Standard Model

### Wide Size

 Different Lids Colors



With Perforated Lid



With Color Options



**Code**      **Outside Dimension**  
**Inside Dimension**  
**L x W x H**

K1-50-110    595 x 365 x 105 mm  
545 x 352 x 100 mm

K1-50-170    595 x 365 x 165 mm  
545 x 352 x 160mm



**Code**      **Outside Dimension**  
**Inside Dimension**  
**L x W x H**

KP1-50-110    595 x 365 x 105 mm  
545 x 352 x 88 mm

KP1-50-170    595 x 365 x 165 mm  
545 x 352 x 148 mm



**Colors**    **Code**

Grey      K1-50-1

Yellow    K1-50-2

Red        K1-50-3

Blue       K1-50-4

Green     K1-50-5

Black      K1-50-6

Purple    K1-50-7

Bronze    K1-50-8


**Dimension**  
585 x 370 mm



### Wire Basket



Mesh: 4x4 mm  
Wire: ø 1 mm



**Code**      **Outside Dimension**  
**L x W x H**

B1-5333-50    533 x 337 x 50 mm

B1-5333-70    533 x 337 x 70 mm


B1-5333-100   533 x 337 x 100 mm

**Code**

B1-5333

**Outside Dimension**  
**L x W x H**

533 x 337 x 5 mm



Wire Basket

### Perforated Basket



Mesh: 4x4 mm  
Wire: ø 1 mm



**Code**      **Outside Dimension**  
**L x W x H**

B2-5333-30    533 x 339 x 40 mm

B2-5333-50    533 x 339 x 60 mm

B2-5333-70    533 x 339 x 80 mm


B2-5333-100   533 x 339 x 110 mm

**Code**

B2-5333

**Outside Dimension**  
**L x W x H**

534 x 343 x 5 mm




Perforated Basket

### Hybrid Basket



Mesh: 4x4 mm  
Sides: ø 3 mm  
Wire: ø 1 mm



**Code**      **Outside Dimension**  
**L x W x H**

B3-5333-50    533 x 337 x 50 mm

B3-5333-70    533 x 337 x 70 mm

B3-5333-100   533 x 337 x 100 mm

**Code**

B3-5333

**Outside Dimension**  
**L x W x H**

533 x 337 x 5 mm



Hybrid Basket

## Standard Model

### Wide Size Baskets

## Standard Model

Flat Size

 Different Lids Colors



With Perforated Lid



With Color Options



**Code**      **Outside Dimension**  
                  **Inside Dimension**  
                  **L x W x H**

|          |  |
|----------|--|
| K1-90-45 | 295 x 270 x 45 mm<br>255 x 255 x 41 mm |
| K1-90-75 | 295 x 270 x 75 mm<br>255 x 255 x 71 mm |



**Colors**      **Code**


|        |         |
|--------|---------|
| Grey   | K1-90-1 |
| Yellow | K1-90-2 |
| Red    | K1-90-3 |
| Blue   | K1-90-4 |
| Green  | K1-90-5 |
| Black  | K1-90-6 |
| Purple | K1-90-7 |
| Bronze | K1-90-8 |

**Dimension**  
295 x 285 mm

### Wire Basket



Mesh: 4x4 mm  
Wire: ø 1 mm



|             |  |
|-------------|--|
| <b>Code</b> | <b>Outside Dimension</b><br><b>L x W x H</b> |
| B1-2525-30  | 255 x 245 x 30 mm                            |
| B1-2525-50  | 255 x 245 x 50 mm                            |



|             |  |
|-------------|--|
| <b>Code</b> | <b>Outside Dimension</b><br><b>L x W x H</b> |
| B1-2525     | 240 x 215 x 5 mm                             |

### Perforated Basket



Mesh: 4x4 mm  
Wire: ø 1 mm



|             |  |
|-------------|--|
| <b>Code</b> | <b>Outside Dimension</b><br><b>L x W x H</b> |
| B2-2525-30  | 242 x 252 x 40 mm                            |
| B2-2525-50  | 242 x 252 x 60 mm                            |



|             |  |
|-------------|--|
| <b>Code</b> | <b>Outside Dimension</b><br><b>L x W x H</b> |
| B2-2525     | 246 x 256 x 5 mm                             |

### Silicone Mat




|             |  |
|-------------|--|
| <b>Code</b> | <b>Outside Dimension</b><br><b>L x W x H</b> |
| A7-2525     | 225 x 235 mm                                 |


## Standard Model

Flat Size Baskets

### Wire Basket



### Perforated Basket



## Standard Model

### Mini Size

 Different Lids Colors



With Perforated Lid



With Color Options



**Code**      **Outside Dimension**  
                  **Inside Dimension**  
                  **L x W x H**


|           |  |
|-----------|--|
| K1-70-40  | 315 x 130 x 40 mm<br>275 x 127 x 30 mm   |
| K1-70-70  | 315 x 130 x 70 mm<br>275 x 127 x 60 mm   |
| K1-70-100 | 315 x 130 x 100 mm<br>275 x 127 x 90 mm  |
| K1-70-130 | 315 x 130 x 130 mm<br>275 x 127 x 120 mm |



| Colors | Code    |
|--------|---------|
| Grey   | K1-70-1 |
| Yellow | K1-70-2 |
| Red    | K1-70-3 |
| Blue   | K1-70-4 |
| Green  | K1-70-5 |
| Black  | K1-70-6 |
| Purple | K1-70-7 |
| Bronze | K1-70-8 |

**Dimension**  
315 x 130 mm

### Wire Basket



Mesh: 4x4 mm  
Wire: ø 1 mm




| Code        | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B1-2712-40  | 270 x 117 x 40 mm              |
| B1-2712-100 | 268 x 117 x 100 mm             |



| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| B1-2712 | 270 x 117 x 5 mm               |

### Perforated Basket



Mesh: 4x4 mm  
Wire: ø 1 mm



| Code        | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B2-2712-40  | 265 x 125 x 50 mm              |
| B2-2712-100 | 265 x 125 x 110 mm             |



| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| B2-2712 | 268 x 128 x 5 mm               |

### Silicone Mat



| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| A7-2712 | 270 x 125 mm                   |


## Standard Model

### Mini Size Baskets

#### Wire Basket



#### Perforated Basket



## Standard Model

Dental Size

 Different Lids Colors



With Perforated Lid



With Color Options



| Code      | Outside Dimension<br>Inside Dimension<br>L x W x H |
|-----------|--|
| K1-60-40  | 320 x 185 x 40 mm<br>280 x 182 x 30 mm             |
| K1-60-70  | 320 x 185 x 70 mm<br>280 x 182 x 40 mm             |
| K1-60-100 | 320 x 185 x 100 mm<br>280 x 182 x 90 mm            |
| K1-60-130 | 320 x 185 x 130 mm<br>280 x 182 x 120 mm           |



| Colors | Code    |
|--------|---------|
| Grey   | K1-60-1 |
| Yellow | K1-60-2 |
| Red    | K1-60-3 |
| Blue   | K1-60-4 |
| Green  | K1-60-5 |
| Black  | K1-60-6 |
| Purple | K1-60-7 |
| Bronze | K1-60-8 |

**Dimension**  
325 x 190 mm

### Wire Basket




Mesh: ø 4x4 mm  
Wire: ø 1 mm



| Code        | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B1-2717-40  | 270 x 170 x 40 mm              |
| B1-2717-100 | 270 x 170 x 100 mm             |

**Code**  
B1-2717  
**Outside Dimension**  
L x W x H  
270 x 170 x 5 mm



Wire Basket

### Perforated Basket



Mesh: ø 4x4 mm  
Wire: ø 1 mm



| Code        | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B2-2717-40  | 275 x 175 x 50 mm              |
| B2-2717-100 | 275 x 175 x 110 mm             |

**Code**  
B2-2717  
**Outside Dimension**  
L x W x H  
278 x 178 x 5 mm



Perforated Basket

## Standard Model

Dental Size Baskets

Wire Basket

Perforated Basket

## Standard Model

1/2 Dental Size

 Different Lids Colors



With Perforated Lid



With Color Options



**Code**      **Outside Dimension**  
                  **Inside Dimension**  
                  **L x W x H**


|           |  |
|-----------|--|
| K1-80-40  | 225 x 160 x 40 mm<br>190 x 157 x 30 mm   |
| K1-80-70  | 225 x 160 x 70 mm<br>190 x 157 x 60 mm   |
| K1-80-100 | 225 x 160 x 100 mm<br>190 x 157 x 90 mm  |
| K1-80-130 | 225 x 160 x 130 mm<br>190 x 157 x 120 mm |



| Colors | Code    |
|--------|---------|
| Grey   | K1-80-1 |
| Yellow | K1-80-2 |
| Red    | K1-80-3 |
| Blue   | K1-80-4 |
| Green  | K1-80-5 |
| Black  | K1-80-6 |
| Purple | K1-80-7 |
| Bronze | K1-80-8 |

**Dimension**  
230 x 175 mm

### Wire Basket



Mesh: ø 4x4 mm  
Wire: ø 1 mm



| Code        | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B1-1714-40  | 175 x 145 x 40 mm              |
| B1-1714-100 | 175 x 145 x 100 mm             |



| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| B1-1714 | 175 x 145 x 5 mm               |

### Perforated Basket



Mesh: ø 4x4 mm  
Wire: ø 1 mm



| Code        | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B2-1714-40  | 175 x 145 x 50 mm              |
| B2-1714-100 | 175 x 145 x 110 mm             |



| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| B2-1714 | 178 x 148 x 5 mm               |

## Standard Model

1/2 Dental Size Baskets

### Wire Basket



### Perforated Basket



Clean. Safe. **Sustainable.**

2 Different Filtering System  
Mechanical Valve System & Barrier System




Different Lids Colors

victorishealth.com

**BARRIER MODEL**

The Barrier Model is an innovative sterilization container designed to ensure maximum safety in CSSD workflows. Its labyrinth structure prevents microorganisms from entering, maintaining sterility for extended periods. This design offers reliable, practical, and scientifically proven protection.



**EASE WITH COLORS  
EXCELLENCE IN WORKFLOW**

Different Lids Colors



**Scientific Foundation**

The principle behind the Barrier mechanism traces back to Louis Pasteur's groundbreaking "Swan Neck Flask" experiment in 1862, which challenged the idea of spontaneous generation. Pasteur designed a flask with an S-shaped neck that allowed air to pass freely but trapped dust and microorganisms in its curves. When the broth inside was boiled and left untouched, no microbial growth occurred, clearly proving that life could not emerge without contamination from existing microorganisms. However, once the flask was tilted and the trapped particles reached the broth, microbial growth appeared almost immediately.

The Barrier Model adopts the same principle through its labyrinth-like system. By forcing air to travel through a complex path, it prevents microorganisms and viruses from penetrating the container. This ensures that sterility is maintained, while still allowing necessary airflow. The result is a reliable and scientifically validated protective system, directly inspired by one of the most important experiments in microbiology.



**REDUCE. REUSE. RENEW.**

Seize Cost Control and Sustainability: Ditch the waste expense and embrace the smart, permanent solution.



**VALVE MODEL**

Victoris Valve Model sterilization containers are designed with an advanced mechanical valve system that removes the need for disposable filters while delivering stable and reliable performance throughout extended use.


Each valve is precisely engineered with integrated gaskets that ensure smooth operation in response to pressure variations during the sterilization process. The secure closure, supported by high-grade silicone gaskets, provides excellent sealing integrity and long-term durability.

Combining innovative design with practicality, Victoris Valve Model containers offer a sustainable and cost-efficient solution for safe and effective sterilization in healthcare environments worldwide.

EASE WITH COLORS EXCELLENCE IN WORKFLOW



**RUNNING PRINCIPLE**



At the beginning of the sterilization cycle, the vacuum phase causes the valves to fully open, allowing the evacuation of dry air from the container.

During the sterilization cycle, the pressure-activated valves automatically open under positive pressure, enabling controlled steam penetration for effective sterilization.

Once the process is complete, and no external forces such as vacuum or pressure are present, the entire system remains sealed, preserving the sterility of its contents.

## Reusable Model

1/1 Full Size

 Different Lids Colors



With Reusable Mechanism



With Color Options



**Non Perforated**

| Code      | Outside Dimension<br>Inside Dimension<br>L x W x H |
|-----------|--|
| K1-10-110 | 595 x 270 x 95 mm<br>544 x 257 x 92 mm             |
| K1-10-135 | 595 x 270 x 115 mm<br>544 x 257 x 113 mm           |
| K1-10-150 | 595 x 270 x 140 mm<br>544 x 257 x 137 mm           |
| K1-10-200 | 595 x 270 x 190 mm<br>544 x 257 x 187 mm           |
| K1-10-260 | 595 x 270 x 250 mm<br>544 x 257 x 247 mm           |



**Valve Lids**

| Colors | Code    |
|--------|---------|
| Grey   | K2-10-1 |
| Yellow | K2-10-2 |
| Red    | K2-10-3 |
| Blue   | K2-10-4 |
| Green  | K2-10-5 |
| Black  | K2-10-6 |
| Purple | K2-10-7 |
| Bronze | K2-10-8 |

**Dimension**  
585 x 285 mm




**Barrier Lids**


| Colors | Code    |
|--------|---------|
| Grey   | K3-10-1 |
| Yellow | K3-10-2 |
| Red    | K3-10-3 |
| Blue   | K3-10-4 |
| Green  | K3-10-5 |
| Black  | K3-10-6 |
| Purple | K3-10-7 |
| Bronze | K3-10-8 |

**Dimension**  
585 x 285 mm

### Wire Basket



Mesh: 4x4 mm  
Wire: ø 1 mm



Box



Lid


| Code        | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B1-5425-30  | 540 x 255 x 30 mm              |
| B1-5425-50  | 540 x 255 x 50 mm              |
| B1-5425-70  | 540 x 255 x 70 mm              |
| B1-5425-100 | 540 x 255 x 100 mm             |

**Code**  
B1-5425  
**Outside Dimension**  
L x W x H  
510 x 240 x 5 mm

### Perforated Basket



Mesh: 4x4 mm  
Wire: ø 1 mm



Box (with 10 mm feet)



Lid


| Code        | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B2-5425-30  | 540 x 253 x 40 mm              |
| B2-5425-50  | 540 x 253 x 60 mm              |
| B2-5425-70  | 540 x 253 x 80 mm              |
| B2-5425-100 | 540 x 253 x 110 mm             |

**Code**  
B2-5425  
**Outside Dimension**  
L x W x H  
540 x 256 x 5 mm

### Hybrid Basket



Mesh: 4x4 mm  
Sides: ø 3 mm  
Wire: ø 1 mm



Box



Lid

| Code        | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B3-5425-50  | 540 x 255 x 50 mm              |
| B3-5425-70  | 540 x 255 x 70 mm              |
| B3-5425-100 | 540 x 255 x 100 mm             |

**Code**  
B3-5425  
**Outside Dimension**  
L x W x H  
540 x 255 x 5 mm

### Silicone Mat



**Code**  
A7-5425  
**Outside Dimension**  
L x W x H  
500 x 235 mm


## Reusable Model

1/1 Full Size Baskets

### Wire Basket



### Perforated Basket



### Hybrid Basket



## Reusable Model

3/4 Quarter Size

 Different Lids Colors



With Reusable Mechanism



With Color Options



**Non Perforated**

| Code      | Outside Dimension<br>Inside Dimension<br>L x W x H |
|-----------|--|
| K1-20-110 | 470 x 270 x 95 mm<br>430 x 260 x 92 mm             |
| K1-20-135 | 470 x 270 x 115 mm<br>430 x 260 x 101 mm           |
| K1-20-150 | 470 x 270 x 140 mm<br>430 x 260 x 137 mm           |
| K1-20-200 | 470 x 270 x 184 mm<br>430 x 260 x 178 mm           |



**Valve Lids**

| Colors | Code    |
|--------|---------|
| Grey   | K2-20-1 |
| Yellow | K2-20-2 |
| Red    | K2-20-3 |
| Blue   | K2-20-4 |
| Green  | K2-20-5 |
| Black  | K2-20-6 |
| Purple | K2-20-7 |
| Bronze | K2-20-8 |

**Dimension**  
465 x 285 mm




**Barrier Lids**


| Colors | Code    |
|--------|---------|
| Grey   | K3-20-1 |
| Yellow | K3-20-2 |
| Red    | K3-20-3 |
| Blue   | K3-20-4 |
| Green  | K3-20-5 |
| Black  | K3-20-6 |
| Purple | K3-20-7 |
| Bronze | K3-20-8 |

**Dimension**  
465 x 285 mm

### Wire Basket



Mesh: 4x4 mm  
Wire: ø 1 mm



Box



Lid

| Code        | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B1-4025-30  | 405 x 255 x 30 mm              |
| B1-4025-50  | 405 x 255 x 50 mm              |
| B1-4025-70  | 405 x 255 x 70 mm              |
| B1-4025-100 | 405 x 255 x 100 mm             |

**Code**  
B1-4025  
**Outside Dimension  
L x W x H**  
375 x 240 x 5 mm

### Perforated Basket



Mesh: 4x4 mm  
Wire: ø 1 mm



Box (with 10 mm feet)



Lid


| Code        | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B2-4025-30  | 405 x 252 x 40 mm              |
| B2-4025-50  | 405 x 252 x 60 mm              |
| B2-4025-70  | 405 x 252 x 80 mm              |
| B2-4025-100 | 405 x 252 x 110 mm             |

**Code**  
B2-4025  
**Outside Dimension  
L x W x H**  
409 x 256 x 5 mm

### Hybrid Basket



Mesh: 4x4 mm  
Sides: ø 3 mm  
Wire: ø 1 mm



Box



Lid

| Code        | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B3-4025-50  | 405 x 245 x 50 mm              |
| B3-4025-70  | 405 x 245 x 70 mm              |
| B3-4025-100 | 405 x 245 x 100 mm             |

**Code**  
B3-4025  
**Outside Dimension  
L x W x H**  
405 x 245 x 5 mm

### Silicone Mat



**Code**  
A7-4025  
**Outside Dimension  
L x W x H**  
370 x 235 mm

## Reusable Model

3/4 Quarter Size Baskets

### Wire Basket



### Perforated Basket



### Hybrid Basket



## Reusable Model

1/2 Half Size

 Different Lids Colors



With Reusable Mechanism



With Color Options



**Non Perforated**

| Code      | Outside Dimension<br>Inside Dimension<br>L x W x H |
|-----------|--|
| K1-30-110 | 300 x 270 x 95 mm<br>255 x 255 x 92 mm             |
| K1-30-135 | 300 x 270 x 115 mm<br>255 x 255 x 113 mm           |
| K1-30-150 | 300 x 270 x 140 mm<br>255 x 255 x 137 mm           |
| K1-30-200 | 300 x 270 x 190 mm<br>255 x 255 x 187 mm           |
| K1-30-260 | 300 x 270 x 250 mm<br>255 x 255 x 247 mm           |



**Valve Lids**

| Colors | Code    |
|--------|---------|
| Grey   | K2-30-1 |
| Yellow | K2-30-2 |
| Red    | K2-30-3 |
| Blue   | K2-30-4 |
| Green  | K2-30-5 |
| Black  | K2-30-6 |
| Purple | K2-30-7 |
| Bronze | K2-30-8 |

**Dimension**  
295 x 285 mm




**Barrier Lids**


| Colors | Code    |
|--------|---------|
| Grey   | K3-30-1 |
| Yellow | K3-30-2 |
| Red    | K3-30-3 |
| Blue   | K3-30-4 |
| Green  | K3-30-5 |
| Black  | K3-30-6 |
| Purple | K3-30-7 |
| Bronze | K3-30-8 |

**Dimension**  
295 x 285 mm

### Wire Basket



Mesh: 4x4 mm  
Wire: ø 1 mm



**Box**



**Lid**

| Code        | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B1-2525-30  | 255 x 245 x 30 mm              |
| B1-2525-50  | 255 x 245 x 50 mm              |
| B1-2525-70  | 255 x 245 x 70 mm              |
| B1-2525-100 | 255 x 245 x 100 mm             |

**Code**  
B1-2525  
**Outside Dimension**  
L x W x H  
240 x 215 x 5 mm


### Perforated Basket



Mesh: 4x4 mm  
Wire: ø 1 mm



**Box (with 10 mm feet)**



**Lid**


| Code        | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B2-2525-30  | 242 x 252 x 40 mm              |
| B2-2525-50  | 242 x 252 x 60 mm              |
| B2-2525-70  | 242 x 252 x 80 mm              |
| B2-2525-100 | 242 x 252 x 110 mm             |

**Code**  
B2-2525  
**Outside Dimension**  
L x W x H  
246 x 256 x 5 mm

### Hybrid Basket



Mesh: 4x4 mm  
Sides: ø 3 mm  
Wire: ø 1 mm



**Box**



**Lid**

| Code        | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B3-2525-50  | 255 x 245 x 50 mm              |
| B3-2525-70  | 255 x 245 x 70 mm              |
| B3-2525-100 | 255 x 245 x 100 mm             |

**Code**  
B3-2525  
**Outside Dimension**  
L x W x H  
245 x 245 x 5 mm

### Silicone Mat



**Code**  
A7-2525  
**Outside Dimension**  
L x W x H  
225 x 235 mm


## Reusable Model

1/2 Half Size Baskets

### Wire Basket



### Perforated Basket



### Hybrid Basket



## Reusable Model

Long Size

 Different Lids Colors



With Reusable Mechanism



With Color Options



**Code**      **Outside Dimension**  
**Inside Dimension**  
**L x W x H**

|           |  |
|-----------|--|
| K1-40-150 | 770 x 270 x 135 mm<br>720 x 260 x 130 mm |
| K1-40-200 | 770 x 270 x 190 mm<br>720 x 260 x 185 mm |



**Colors**    **Code**

|        |         |
|--------|---------|
| Grey   | K2-40-1 |
| Yellow | K2-40-2 |
| Red    | K2-40-3 |
| Blue   | K2-40-4 |
| Green  | K2-40-5 |
| Black  | K2-40-6 |
| Purple | K2-40-7 |
| Bronze | K2-40-8 |

**Dimension**  
770 x 270 mm





**Colors**    **Code**

|        |         |
|--------|---------|
| Grey   | K3-40-1 |
| Yellow | K3-40-2 |
| Red    | K3-40-3 |
| Blue   | K3-40-4 |
| Green  | K3-40-5 |
| Black  | K3-40-6 |
| Purple | K3-40-7 |
| Bronze | K3-40-8 |

**Dimension**  
770 x 270 mm

### Wire Basket




|             |  |
|-------------|--|
| <b>Code</b> | <b>Outside Dimension</b><br><b>L x W x H</b> |
| B1-6625-50  | 665 x 250 x 50 mm                            |
| B1-6625-70  | 665 x 250 x 70 mm                            |
| B1-6625-100 | 665 x 250 x 100 mm                           |

Mesh: 4x4 mm  
Wire: ø 1 mm



|             |  |
|-------------|--|
| <b>Code</b> | <b>Outside Dimension</b><br><b>L x W x H</b> |
| B1-6625     | 635 x 235 x 5 mm                             |

### Perforated Basket




|             |  |
|-------------|--|
| <b>Code</b> | <b>Outside Dimension</b><br><b>L x W x H</b> |
| B2-6625-30  | 665 x 253 x 40 mm                            |
| B2-6625-50  | 665 x 253 x 60 mm                            |
| B2-6625-70  | 665 x 253 x 80 mm                            |
| B2-6625-100 | 665 x 253 x 110 mm                           |

Mesh: 4x4 mm  
Wire: ø 1 mm



|             |  |
|-------------|--|
| <b>Code</b> | <b>Outside Dimension</b><br><b>L x W x H</b> |
| B2-6625     | 668 x 256 x 5 mm                             |

### Hybrid Basket



|             |  |
|-------------|--|
| <b>Code</b> | <b>Outside Dimension</b><br><b>L x W x H</b> |
| B3-6625-50  | 665 x 245 x 50 mm                            |
| B3-6625-70  | 665 x 245 x 70 mm                            |
| B3-6625-100 | 665 x 245 x 100 mm                           |

Mesh: 4x4 mm  
Sides: ø 3 mm  
Wire: ø 1 mm



|             |  |
|-------------|--|
| <b>Code</b> | <b>Outside Dimension</b><br><b>L x W x H</b> |
| B3-6625     | 668 x 256 x 5 mm                             |


## Reusable Model

Long Size Baskets

### Wire Basket



### Perforated Basket



### Hybrid Basket



## Reusable Model

Wide Size


 Different Lids Colors



With Reusable Mechanism



With Color Options



| Code      | Outside Dimension<br>Inside Dimension<br>L x W x H |
|-----------|--|
| K1-50-110 | 595 x 365 x 105 mm<br>545 x 352 x 100 mm           |
| K1-50-170 | 595 x 365 x 165 mm<br>545 x 352 x 160 mm           |



| Colors | Code    |
|--------|---------|
| Grey   | K2-50-1 |
| Yellow | K2-50-2 |
| Red    | K2-50-3 |
| Blue   | K2-50-4 |
| Green  | K2-50-5 |
| Black  | K2-50-6 |
| Purple | K2-50-7 |
| Bronze | K2-50-8 |



**Dimension**  
585 x 370 mm



| Colors | Code    |
|--------|---------|
| Grey   | K3-50-1 |
| Yellow | K3-50-2 |
| Red    | K3-50-3 |
| Blue   | K3-50-4 |
| Green  | K3-50-5 |
| Black  | K3-50-6 |
| Purple | K3-50-7 |
| Bronze | K3-50-8 |

**Dimension**  
585 x 370 mm

### Wire Basket



| Code        | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B1-5333-50  | 533 x 337 x 50 mm              |
| B1-5333-70  | 533 x 337 x 70 mm              |
| B1-5333-100 | 533 x 337 x 100 mm             |

Mesh: 4x4 mm  
Wire: ø 1 mm



| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| B1-5333 | 533 x 337 x 5 mm               |

### Perforated Basket




| Code        | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B2-5333-30  | 533 x 339 x 40 mm              |
| B2-5333-50  | 533 x 339 x 60 mm              |
| B2-5333-70  | 533 x 339 x 80 mm              |
| B2-5333-100 | 533 x 339 x 110 mm             |

Mesh: 4x4 mm  
Wire: ø 1 mm



| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| B2-5333 | 534 x 343 x 5 mm               |

### Hybrid Basket



| Code        | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B3-5333-50  | 533 x 337 x 50 mm              |
| B3-5333-70  | 533 x 337 x 70 mm              |
| B3-5333-100 | 533 x 337 x 100 mm             |

Mesh: 4x4 mm  
Sides: ø 3 mm  
Wire: ø 1 mm




| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| B3-5333 | 533 x 337 x 5 mm               |

## Reusable Model

Wide Size Baskets

### Wire Basket



### Perforated Basket



### Hybrid Basket



## Reusable Model

Flat Size

 Different Lids Colors



With Reusable Mechanism



With Color Options



**Code**      **Outside Dimension**  
                  **Inside Dimension**  
                  **L x W x H**

|          |  |
|----------|--|
| K1-90-75 | 295 x 270 x 75 mm<br>255 x 255 x 71 mm |
|----------|--|



**Colors**    **Code**

|        |         |
|--------|---------|
| Grey   | K3-90-1 |
| Yellow | K3-90-2 |
| Red    | K3-90-3 |
| Blue   | K3-90-4 |
| Green  | K3-90-5 |
| Black  | K3-90-6 |
| Purple | K3-90-7 |
| Bronze | K3-90-8 |

**Dimension**  
295 x 285 mm

### Wire Basket



Mesh: 4x4 mm  
Wire: ø 1 mm




**Code**      **Outside Dimension**  
                  **L x W x H**

|            |                   |
|------------|-------------------|
| B1-2525-30 | 255 x 245 x 30 mm |
| B1-2525-50 | 255 x 245 x 50 mm |



**Code**  
B1-2525  
**Outside Dimension**  
**L x W x H**  
240 x 215 x 5 mm

### Perforated Basket



Mesh: 4x4 mm  
Wire: ø 1 mm



**Code**      **Outside Dimension**  
                  **L x W x H**

|            |                   |
|------------|-------------------|
| B2-2525-30 | 242 x 252 x 40 mm |
| B2-2525-50 | 242 x 252 x 60 mm |



**Code**  
B2-2525  
**Outside Dimension**  
**L x W x H**  
246 x 256 x 5 mm

### Silicone Mat




**Code**  
A7-2525  
**Outside Dimension**  
**L x W x H**  
225 x 235 mm


## Reusable Model

Flat Size Baskets

### Wire Basket



### Perforated Basket



## Reusable Model

### Mini Size

 Different Lids Colors



With Reusable Mechanism



With Color Options



**Code**      **Outside Dimension**  
                  **Inside Dimension**  
                  **L x W x H**


|           |  |
|-----------|--|
| K1-70-70  | 315 x 130 x 70 mm<br>275 x 127 x 60 mm   |
| K1-70-100 | 315 x 130 x 100 mm<br>275 x 127 x 90 mm  |
| K1-70-130 | 315 x 130 x 130 mm<br>275 x 127 x 120 mm |



| Colors | Code    |
|--------|---------|
| Grey   | K3-70-1 |
| Yellow | K3-70-2 |
| Red    | K3-70-3 |
| Blue   | K3-70-4 |
| Green  | K3-70-5 |
| Black  | K3-70-6 |
| Purple | K3-70-7 |
| Bronze | K3-70-8 |

**Dimension**  
315 x 130 mm

### Wire Basket



Mesh: 4x4 mm  
Wire: ø 1 mm



| Code        | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B1-2712-40  | 270 x 117 x 40 mm              |
| B1-2712-100 | 268 x 117 x 100 mm             |



| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| B1-2712 | 270 x 117 x 5 mm               |

### Perforated Basket



Mesh: 4x4 mm  
Wire: ø 1 mm



| Code        | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B2-2712-40  | 265 x 125 x 50 mm              |
| B2-2712-100 | 265 x 125 x 110 mm             |



| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| B2-2712 | 268 x 128 x 5 mm               |

### Silicone Mat



| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| A7-2712 | 270 x 125 mm                   |


## Reusable Model

### Mini Size Baskets

#### Wire Basket



#### Perforated Basket



## Reusable Model

Dental Size

 Different Lids Colors



With Reusable Mechanism



With Color Options



**Code**      **Outside Dimension**  
                  **Inside Dimension**  
                  **L x W x H**

|           |  |
|-----------|--|
| K1-60-70  | 320 x 185 x 70 mm<br>280 x 182 x 40 mm   |
| K1-60-100 | 320 x 185 x 100 mm<br>280 x 182 x 90 mm  |
| K1-60-130 | 320 x 185 x 130 mm<br>280 x 182 x 120 mm |



**Colors**      **Code**


|        |         |
|--------|---------|
| Grey   | K3-60-1 |
| Yellow | K3-60-2 |
| Red    | K3-60-3 |
| Blue   | K3-60-4 |
| Green  | K3-60-5 |
| Black  | K3-60-6 |
| Purple | K3-60-7 |
| Bronze | K3-60-8 |

**Dimension**  
 320 x 190 mm

### Wire Basket



Mesh: ø 4x4 mm  
 Wire: ø 1 mm



**Code**      **Outside Dimension**  
                  **L x W x H**

|             |                    |
|-------------|--------------------|
| B1-2717-40  | 270 x 170 x 40 mm  |
| B1-2717-100 | 270 x 170 x 100 mm |



**Code**  
 B1-2717  
**Outside Dimension**  
**L x W x H**  
 270 x 170 x 5 mm

### Perforated Basket



Mesh: ø 4x4 mm  
 Wire: ø 1 mm



**Code**      **Outside Dimension**  
                  **L x W x H**

|             |                    |
|-------------|--------------------|
| B2-2717-40  | 275 x 175 x 50 mm  |
| B2-2717-100 | 275 x 175 x 110 mm |




**Code**  
 B2-2717  
**Outside Dimension**  
**L x W x H**  
 278 x 178 x 5 mm

## Reusable Model

Dental Size Baskets

### Wire Basket



### Perforated Basket



## Reusable Model

1/2 Dental Size

 Different Lids Colors



With Reusable Mechanism



With Color Options



Non Perforated

**Code**      **Outside Dimension**  
**Inside Dimension**  
**L x W x H**

|           |  |
|-----------|--|
| K1-80-70  | 225 x 160 x 70 mm<br>190 x 157 x 60 mm   |
| K1-80-100 | 225 x 160 x 100 mm<br>190 x 157 x 90 mm  |
| K1-80-130 | 225 x 160 x 130 mm<br>190 x 157 x 120 mm |




Lids

| Colors | Code    |
|--------|---------|
| Grey   | K3-80-1 |
| Yellow | K3-80-2 |
| Red    | K3-80-3 |
| Blue   | K3-80-4 |
| Green  | K3-80-5 |
| Black  | K3-80-6 |
| Purple | K3-80-7 |
| Bronze | K3-80-8 |

**Dimension**  
230 x 175 mm

### Wire Basket



Mesh: ø 4x4 mm  
Wire: ø 1 mm



Box

**Code**      **Outside Dimension**  
**L x W x H**

|             |                    |
|-------------|--------------------|
| B1-1714-40  | 175 x 145 x 40 mm  |
| B1-1714-100 | 175 x 145 x 100 mm |



Lid

**Code**  
B1-1714  
**Outside Dimension**  
**L x W x H**  
175 x 145 x 5 mm

### Perforated Basket



Mesh: ø 4x4 mm  
Wire: ø 1 mm



Box (with 10 mm feet)

**Code**      **Outside Dimension**  
**L x W x H**

|             |                    |
|-------------|--------------------|
| B2-1714-40  | 175 x 145 x 50 mm  |
| B2-1714-100 | 175 x 145 x 110 mm |



Lid

**Code**  
B2-1714  
**Outside Dimension**  
**L x W x H**  
178 x 148 x 5 mm

## Reusable Model

1/2 Dental Size Baskets

### Wire Basket



### Perforated Basket



Organize.

Sterilize.

Excel.

Holding the Standard




Three Different Types of Basket Models by Construction/Material




Perforated Sheet  
Ø 1 mm  
□ 4 mm



Wire Mesh  
Ø 1 mm  
□ 4 mm



Hybrith  
Ø 1 mm  
□ 4 mm  
Side Mesh  
Ø 3 mm



**Half Size**

| Box Code    | Outside Dimension L x W x H |
|-------------|-----------------------------|
| B1-2525-30  | 255 x 245 x 30 mm           |
| B1-2525-50  | 255 x 245 x 50 mm           |
| B1-2525-70  | 255 x 245 x 70 mm           |
| B1-2525-100 | 255 x 245 x 100 mm          |

| Lid Code | Outside Dimension L x W x H |
|----------|-----------------------------|
| B1-2525  | 240 x 215 x 5 mm            |

**Long Size**

| Box Code    | Outside Dimension L x W x H |
|-------------|-----------------------------|
| B1-6625-50  | 665 x 250 x 50 mm           |
| B1-6625-70  | 665 x 250 x 70 mm           |
| B1-6625-100 | 665 x 250 x 100 mm          |
| -           | -                           |

| Lid Code | Outside Dimension L x W x H |
|----------|-----------------------------|
| B1-6625  | 635 x 235 x 5 mm            |

**Wide Size**

| Box Code    | Outside Dimension L x W x H |
|-------------|-----------------------------|
| B1-5333-50  | 533 x 337 x 50 mm           |
| B1-5333-70  | 533 x 337 x 70 mm           |
| B1-5333-100 | 533 x 337 x 100 mm          |
| -           | -                           |

| Lid Code | Outside Dimension L x W x H |
|----------|-----------------------------|
| B1-5333  | 533 x 337 x 5 mm            |

**Full Size**

| Box Code    | Outside Dimension L x W x H |
|-------------|-----------------------------|
| B1-5425-30  | 540 x 255 x 30 mm           |
| B1-5425-50  | 540 x 255 x 50 mm           |
| B1-5425-70  | 540 x 255 x 70 mm           |
| B1-5425-100 | 540 x 255 x 100 mm          |

| Lid Code | Outside Dimension L x W x H |
|----------|-----------------------------|
| B1-5425  | 510 x 240 x 5 mm            |

**Din Size**

| Box Code    | Outside Dimension L x W x H |
|-------------|-----------------------------|
| B1-4825-30  | 485 x 255 x 30 mm           |
| B1-4825-50  | 485 x 255 x 50 mm           |
| B1-4825-70  | 485 x 255 x 70 mm           |
| B1-4825-100 | 485 x 255 x 100 mm          |

| Lid Code | Outside Dimension L x W x H |
|----------|-----------------------------|
| B1-4825  | 455 x 240 x 5 mm            |

**Quarter Size**

| Box Code    | Outside Dimension L x W x H |
|-------------|-----------------------------|
| B1-4025-30  | 405 x 255 x 30 mm           |
| B1-4025-50  | 405 x 255 x 50 mm           |
| B1-4025-70  | 405 x 255 x 70 mm           |
| B1-4025-100 | 405 x 255 x 100 mm          |

| Lid Code | Outside Dimension L x W x H |
|----------|-----------------------------|
| B1-4025  | 375 x 240 x 5 mm            |

**Dental Size**

| Box Code    | Outside Dimension L x W x H |
|-------------|-----------------------------|
| B1-2717-40  | 270 x 170 x 40 mm           |
| B1-2717-100 | 270 x 170 x 100 mm          |
| -           | -                           |
| -           | -                           |

| Lid Code | Outside Dimension L x W x H |
|----------|-----------------------------|
| B1-2717  | 270 x 170 x 5 mm            |

**Mini Size**


| Box Code    | Outside Dimension L x W x H |
|-------------|-----------------------------|
| B1-2712-40  | 270 x 117 x 40 mm           |
| B1-2712-100 | 270 x 117 x 100 mm          |
| -           | -                           |
| -           | -                           |

| Lid Code | Outside Dimension L x W x H |
|----------|-----------------------------|
| B1-2712  | 270 x 117 x 5 mm            |

**Dental 1/2**

| Box Code    | Outside Dimension L x W x H |
|-------------|-----------------------------|
| B1-1714-40  | 175 x 145 x 40 mm           |
| B1-1714-100 | 175 x 145 x 100 mm          |
| -           | -                           |
| -           | -                           |

| Lid Code | Outside Dimension L x W x H |
|----------|-----------------------------|
| B1-1714  | 175 x 145 x 5 mm            |



Half Size



| Box Code    | Outside Dimension L x W x H |
|-------------|-----------------------------|
| B2-2525-30  | 242 x 252 x 40 mm           |
| B2-2525-50  | 242 x 252 x 60 mm           |
| B2-2525-70  | 242 x 252 x 80 mm           |
| B2-2525-100 | 242 x 252 x 110 mm          |

| Lid Code | Outside Dimension L x W x H |
|----------|-----------------------------|
| B2-2525  | 246 x 256 x 5 mm            |

Long Size



| Box Code    | Outside Dimension L x W x H |
|-------------|-----------------------------|
| B2-6625-30  | 665 x 253 x 40 mm           |
| B2-6625-50  | 665 x 253 x 60 mm           |
| B2-6625-70  | 665 x 253 x 80 mm           |
| B2-6625-100 | 665 x 253 x 110 mm          |

| Lid Code | Outside Dimension L x W x H |
|----------|-----------------------------|
| B2-6625  | 668 x 256 x 5 mm            |

Wide Size



| Box Code    | Outside Dimension L x W x H |
|-------------|-----------------------------|
| B2-5333-30  | 533 x 339 x 40 mm           |
| B2-5333-50  | 533 x 339 x 60 mm           |
| B2-5333-70  | 533 x 339 x 80 mm           |
| B2-5333-100 | 533 x 339 x 110 mm          |

| Lid Code | Outside Dimension L x W x H |
|----------|-----------------------------|
| B2-5333  | 536 x 343 x 5 mm            |

Full Size



| Box Code    | Outside Dimension L x W x H |
|-------------|-----------------------------|
| B2-5425-30  | 540 x 255 x 40 mm           |
| B2-5425-50  | 540 x 255 x 60 mm           |
| B2-5425-70  | 540 x 255 x 80 mm           |
| B2-5425-100 | 540 x 255 x 110 mm          |

| Lid Code | Outside Dimension L x W x H |
|----------|-----------------------------|
| B2-5425  | 543 x 256 x 5 mm            |

Din Size



| Box Code    | Outside Dimension L x W x H |
|-------------|-----------------------------|
| B2-4825-30  | 484 x 252 x 40 mm           |
| B2-4825-50  | 484 x 252 x 60 mm           |
| B2-4825-70  | 484 x 252 x 80 mm           |
| B2-4825-100 | 484 x 252 x 110 mm          |

| Lid Code | Outside Dimension L x W x H |
|----------|-----------------------------|
| B2-4825  | 488 x 256 x 5 mm            |

Quarter Size



| Box Code    | Outside Dimension L x W x H |
|-------------|-----------------------------|
| B2-4025-30  | 405 x 252 x 40 mm           |
| B2-4025-50  | 405 x 252 x 60 mm           |
| B2-4025-70  | 405 x 252 x 80 mm           |
| B2-4025-100 | 405 x 252 x 110 mm          |

| Lid Code | Outside Dimension L x W x H |
|----------|-----------------------------|
| B2-4025  | 409 x 256 x 5 mm            |

Dental Size



| Box Code    | Outside Dimension L x W x H |
|-------------|-----------------------------|
| B2-2717-40  | 275 x 175 x 50 mm           |
| B2-2717-100 | 275 x 175 x 110 mm          |
| -           | -                           |
| -           | -                           |

| Lid Code | Outside Dimension L x W x H |
|----------|-----------------------------|
| B2-2717  | 278 x 178 x 5 mm            |

Mini Size



| Box Code    | Outside Dimension L x W x H |
|-------------|-----------------------------|
| B2-2712-40  | 265 x 125 x 50 mm           |
| B2-2712-100 | 265 x 125 x 110 mm          |
| -           | -                           |
| -           | -                           |


| Lid Code | Outside Dimension L x W x H |
|----------|-----------------------------|
| B2-2712  | 268 x 128 x 5 mm            |

Dental 1/2



| Box Code    | Outside Dimension L x W x H |
|-------------|-----------------------------|
| B2-1714-40  | 175 x 145 x 50 mm           |
| B2-1714-100 | 175 x 145 x 110 mm          |
| -           | -                           |
| -           | -                           |



| Lid Code | Outside Dimension L x W x H |
|----------|-----------------------------|
| B2-1714  | 178 x 148 x 5 mm            |




| Half Size       |                                    | Long Size       |                                    | Wide Size       |                                    |
|-----------------|------------------------------------|-----------------|------------------------------------|-----------------|------------------------------------|
|                 |                                    |                 |                                    |                 |                                    |
| <b>Box Code</b> | <b>Outside Dimension L x W x H</b> | <b>Box Code</b> | <b>Outside Dimension L x W x H</b> | <b>Box Code</b> | <b>Outside Dimension L x W x H</b> |
| B3-2525-50      | 245 x 245 x 50 mm                  | B3-6625-50      | 665 x 245 x 50 mm                  | B3-5333-50      | 533 x 337 x 50 mm                  |
| B3-2525-70      | 245 x 245 x 70 mm                  | B3-6625-70      | 665 x 245 x 70 mm                  | B3-5333-70      | 533 x 337 x 70 mm                  |
| B3-2525-100     | 245 x 245 x 100 mm                 | B3-6625-100     | 665 x 245 x 100 mm                 | B3-5333-100     | 533 x 337 x 100 mm                 |
| <b>Lid Code</b> | <b>Outside Dimension L x W x H</b> | <b>Lid Code</b> | <b>Outside Dimension L x W x H</b> | <b>Lid Code</b> | <b>Outside Dimension L x W x H</b> |
| B3-2525         | 245 x 245 x 5 mm                   | B3-6625         | 665 x 245 x 5 mm                   | B3-5333         | 533 x 337 x 5 mm                   |

Fine Mesh Baskets

| Full Size       |                                    | Din Size        |                                    | Quarter Size    |                                    |
|-----------------|------------------------------------|-----------------|------------------------------------|-----------------|------------------------------------|
|                 |                                    |                 |                                    |                 |                                    |
| <b>Box Code</b> | <b>Outside Dimension L x W x H</b> | <b>Box Code</b> | <b>Outside Dimension L x W x H</b> | <b>Box Code</b> | <b>Outside Dimension L x W x H</b> |
| B3-5425-50      | 540 x 255 x 50 mm                  | B3-4825-50      | 484 x 252 x 40 mm                  | B3-4025-50      | 405 x 252 x 45 mm                  |
| B3-5425-70      | 540 x 255 x 70 mm                  | B3-4825-70      | 484 x 252 x 70 mm                  | B3-4025-70      | 405 x 252 x 70 mm                  |
| B3-5425-100     | 540 x 255 x 100 mm                 | B3-4825-100     | 484 x 252 x 100 mm                 | B3-4025-100     | 405 x 252 x 100 mm                 |
| <b>Lid Code</b> | <b>Outside Dimension L x W x H</b> | <b>Lid Code</b> | <b>Outside Dimension L x W x H</b> | <b>Lid Code</b> | <b>Outside Dimension L x W x H</b> |
| B3-5425         | 540 x 245 x 5 mm                   | B3-4825         | 485 x 245 x 5 mm                   | B3-4025         | 405 x 245 x 5 mm                   |



| <b>Box Code</b> | <b>Outside Dimension L x W x H</b> | <b>Box Code</b> | <b>Outside Dimension L x W x H</b> |
|-----------------|------------------------------------|-----------------|------------------------------------|
| B4-33-30        | 30 x 30 x 30 mm                    | B4-88-40        | 80 x 80 x 40 mm                    |
| B4-44-20        | 40 x 40 x 20 mm                    | B4-88-60        | 80 x 80 x 60 mm                    |
| B4-66-20        | 60 x 60 x 20 mm                    | B4-1212-30      | 120 x 120 x 30 mm                  |
| B4-66-40        | 60 x 60 x 40 mm                    | B4-1212-60      | 120 x 120 x 60 mm                  |
| B4-84-40        | 80 x 40 x 40 mm                    | B4-2424-60      | 240 x 240 x 60 mm                  |



 Bottom With Lid  
W/O Feet  
Wire:  $\varnothing$  1 mm  
Mesh: 4x4 mm

| Code       | Dimension<br>L x W x H | Code      | Dimension<br>L x W x H |
|------------|------------------------|-----------|------------------------|
| B2-2718-32 | 275 x 180 x 32 mm      | B2-188-32 | 188 x 86 x 32 mm       |
| B2-2718-42 | 275 x 180 x 42 mm      | B2-144-28 | 146 x 46 x 28 mm       |
| B2-247-30  | 240 x 76 x 30 mm       | B2-84-28  | 86 x 46 x 28 mm        |
| B2-1913-34 | 190 x 130 x 34 mm      |           |                        |

STU - SPRI Baskets

Wire Basket  
Single Scope

**Code**  
B5-298-55-S

**Dimension  
L x W x H**  
290 x 80 x 55 mm

**Code**  
B5-468-55-S

**Dimension  
L x W x H**  
460 x 80 x 55 mm

Wire Basket  
Double Scope

**Code**  
B5-298-55-D

**Dimension  
L x W x H**  
290 x 80 x 55 mm

**Code**  
B5-468-55-D

**Dimension  
L x W x H**  
460 x 80 x 55 mm

Perforated Basket  
Single Scope

**Code**  
B5-278-44-S

**Dimension  
L x W x H**  
290 x 80 x 55 mm

**Code**  
B5-458-44-S

**Dimension  
L x W x H**  
452 x 83 x 44 mm

Perforated Basket  
Double Scope

**Code**  
B5-2710-44-D

**Dimension  
L x W x H**  
290 x 80 x 55 mm


**Code**  
B5-4510-44-D

**Dimension  
L x W x H**  
452 x 83 x 44 mm



STU Baskets

| Box Code    | Outside Dimension<br>L x W x H | Details   |
|-------------|--------------------------------|---|
| B7-5728-135 | 575 x 280 x 135 mm             | Upper Frame: $\varnothing$ 6 mm                                     |
| B7-5728-265 | 575 x 280 x 265 mm             | Body Wires: $\varnothing$ 2 & 3 mm<br>Base Frame $\varnothing$ 5 mm |



SPRI Baskets

| Box Code    | Outside Dimension<br>L x W x H | Details                        |
|-------------|--------------------------------|--------------------------------|
| B8-6030-40  | 600x300x40 mm                  | Main Frame: $\varnothing$ 6 mm |
| B8-6030-100 | 600x300x100 mm                 | Body Wire: $\varnothing$ 4 mm  |
| B8-6030-150 | 600x300x150 mm                 | Base Frame $\varnothing$ 5 mm  |
| B8-6040-100 | 600x400x100 mm                 | Single Frame Stackable         |
| B8-6040-150 | 600x400x150 mm                 |                                |
| B8-6040-200 | 600x400x200 mm                 |                                |



| Paper Filters   |  | PTFE Teflon Filters   |  |
|---|--|---|--|
|   |  |   |  |
| <b>Code</b> A1-10-1<br><b>For:</b> 1/1, 3/4, 1/2, Long, Wide<br><b>Size:</b> Ø 185 mm | <b>Code</b> A1-11-1<br><b>For:</b> Small Size Containers<br><b>Size:</b> Ø 88 mm | <b>Code</b> A1-20-1<br><b>For:</b> 1/1, 3/4, 1/2, Long, Wide<br><b>Size:</b> Ø 190 mm | <b>Code</b> A1-21-1<br><b>For:</b> Small Size Containers<br><b>Size:</b> Ø 88 mm |
| <b>Material</b> 73g/m <sup>2</sup> Medical Grade Paper                                |  | <b>Material</b> 1 mm (± 0,25) PTFE  |  |
| <b>Package</b> 100 pcs  |  | <b>Package</b> 1 Piece  |  |
| <b>Sterilization</b> Steam  |  | <b>Sterilization</b> Steam & EtO  |  |
| <b>Cycle</b> Single Use   |  | <b>Cycle</b> 2000 Cycle Reusable  |  |

| Polymer Security Seals  |  | Polymer Security Seals With Indicator   |  |
|---|--|---|--|
|   |  |   |  |
| <b>Code</b> A2-10-1<br><b>Package</b> 100 pcs<br><b>Sterilization</b> Steam<br><b>Cycle</b> Single Use<br><b>Color</b> Blue | <b>Code</b> A2-20-1<br><b>Package</b> 100 pcs<br><b>Sterilization</b> Steam<br><b>Cycle</b> Single Use<br><b>Color</b> Red | <b>Code</b> A2-10-2<br><b>Package</b> 100 pcs<br><b>Sterilization</b> Steam<br><b>Cycle</b> Single Use<br><b>Color</b> Blue | <b>Code</b> A2-20-2<br><b>Package</b> 100 pcs<br><b>Sterilization</b> Steam<br><b>Cycle</b> Single Use<br><b>Color</b> Red |



| Identification Labels             |                     |                            |                     | Filter Retainers   |  |
|-----------------------------------|---------------------|----------------------------|---------------------|--|--|
|                                   | <b>Code</b> A5-10-1 |                            | <b>Code</b> A5-10-5 |  |  |
|                                   | <b>Code</b> A5-10-2 |                            | <b>Code</b> A5-10-6 |  |  |
|                                   | <b>Code</b> A5-10-3 |                            | <b>Code</b> A5-10-7 |  |  |
|                                   | <b>Code</b> A5-10-4 |                            | <b>Code</b> A5-10-8 |  |  |
| <b>Material</b> 1,5 mm Aluminum   |                     | <b>Package</b> 1 Piece     |                     | <b>Code</b> A6-10-1<br><b>For:</b> 1/1, 3/4, 1/2, Long, Wide | <b>Code</b> A6-10-2<br><b>For:</b> Small Size Containers |
| <b>Sterilization</b> Steam        |                     | <b>Sterilization</b> Steam |                     | <b>Material</b> Stainless Steel                              |  |
| <b>Marking</b> Maximum 24 Letters |                     |                            |                     | <b>Package</b> 1 piece                                       |  |

| Paper Labels               |                                       |                            |                                       |
|----------------------------|---------------------------------------|----------------------------|---------------------------------------|
|                            | <b>Code</b> A3-10-1<br>with indicator |                            | <b>Code</b> A3-10-2<br>with indicator |
| <b>Material</b> Paper      |                                       | <b>Material</b> Paper      |                                       |
| <b>Size</b> Big & Small    |                                       | <b>Size</b> Big & Small    |                                       |
| <b>Package</b> 100 pcs     |                                       | <b>Package</b> 100 pcs     |                                       |
| <b>Sterilization</b> Steam |                                       | <b>Sterilization</b> Steam |                                       |

| Silicone Mats & Nets |                     |  |  |
|----------------------|---------------------|--|--|
|                      | <b>Code</b> A7-5425 | <b>Size</b> 500 X 235 mm               |  |
|                      | <b>Code</b> A7-4025 | <b>Size</b> 370 X 235 mm               |  |
|                      | <b>Code</b> A7-2525 | <b>Size</b> 225 x 235 mm               |  |
|                      | <b>Code</b> A7-2712 | <b>Size</b> 270 x 125 mm               |  |
|                      | <b>Code</b> A4-5425 | <b>Size</b> 520 x 242 mm               |  |
|                      | <b>Code</b> A4-4025 | <b>Size</b> 380 x 242 mm               |  |
|                      | <b>Code</b> A4-2525 | <b>Size</b> 255 x 250 mm               |  |
|                      |                     | <b>Material</b> Medical Grade Silicone |  |
|                      |                     | <b>Package</b> 1 Piece                 |  |
|                      |                     | <b>Sterilization</b> Steam             |  |
|                      |                     | <b>Color</b> Light Blue                |  |

| Silicone Gasket |                      |                         |  |                    |                      |                        |
|-----------------|----------------------|-------------------------|--|--------------------|----------------------|------------------------|
|                 | <b>Code</b> A11-10-1 | <b>Size</b> 1/1 Full    | <b>Code</b> A11-40-1                           | <b>Size</b> Long   | <b>Code</b> A11-70-1 | <b>Size</b> Mini       |
|                 | <b>Code</b> A11-20-1 | <b>Size</b> 3/4 Quarter | <b>Code</b> A11-50-1                           | <b>Size</b> Wide   | <b>Code</b> A11-80-1 | <b>Size</b> 1/2 Dental |
|                 | <b>Code</b> A11-30-1 | <b>Size</b> 1/2 Half    | <b>Code</b> A11-60-1                           | <b>Size</b> Dental | <b>Code</b> A11-90-1 | <b>Size</b> Flat       |
|                 |                      |                         | <b>Material</b> Medical Grade Sealing Silicone |                    |                      |                        |
|                 |                      |                         | <b>Package</b> 1 Piece                         |                    |                      |                        |
|                 |                      |                         | <b>Shape</b> D shape, Jointless                |                    |                      |                        |
|                 |                      |                         | <b>Feature</b> Place & Go, No Glue Required    |                    |                      |                        |



**Code T1-10-1**  
 Material : Stainless Steel  
 Feet : 4 wheels - 2 braked  
 Package : 1 Piece  
 Dimensions: 600 x 860 x 900 mm (W x L x H)



**Code T1-10-2**  
 Material : Stainless Steel  
 Feet : 4 wheels - 2 braked  
 Package : 1 Piece  
 Dimensions: 550 x 860 x 900 mm (W x L x H)



**Code T1-10-3**  
 Material : Stainless Steel  
 Feet : 4 wheels - 2 braked  
 Package : 1 Piece  
 Dimensions: 600 x 850 x 900 mm (W x L x H)



**Code T1-10-4**  
 Material : Stainless Steel  
 Feet : 4 wheels - 2 braked  
 Package : 1 Piece  
 Dimensions: 550 x 860 x 900 mm (W x L x H)



**Code T1-10-5**  
 Material : Stainless Steel  
 Feet : 4 Adjustable Leveling Feet  
 Package : 1 Piece  
 Dimensions: 600 x 1400 x 850 mm (W x L x H)



**Code T1-10-6**  
 Material : Stainless Steel  
 Feet : 4 wheels - 2 braked  
 Package : 1 Piece  
 Dimensions: 700 x 360 x 1350 mm (W x L x H)



**Code T1-10-7**  
 Material : Stainless Steel  
 Feet : 4 wheels - 2 braked  
 Package : 1 Piece  
 Dimensions: 700 x 1080 x 1350 mm (W x L x H)



**Code T1-10-8**  
 Material : Stainless Steel  
 Feet : 4 wheels - 2 braked  
 Package : 1 Piece - including 3 baskets  
 Dimensions: 400 x 600 x 1700 (W x L x H)

| Code        | Page                   | Code        | Page               | Code         | Page       | Code      | Page   |
|-------------|------------------------|-------------|--------------------|--------------|------------|-----------|--------|
| A1-10-1     | 68                     | B1-6625-100 | 25, 47, 61         | B3-5425-50   | 19, 41, 65 | K1-30-6   | 22     |
| A11-10-1    | 69                     | B1-6625-50  | 25, 47, 61         | B3-5425-70   | 19, 41, 65 | K1-30-7   | 22     |
| A1-11-1     | 68                     | B1-6625-70  | 25, 47, 61         | B3-6625      | 25, 47, 65 | K1-30-8   | 22     |
| A11-20-1    | 69                     | B2-144-28   | 67                 | B3-6625-100  | 25, 47, 65 | K1-40-1   | 24     |
| A11-30-1    | 69                     | B2-1714     | 35, 57, 63         | B3-6625-50   | 25, 47, 65 | K1-40-150 | 24, 46 |
| A11-40-1    | 69                     | B2-1714-100 | 35, 57, 63         | B3-6625-70   | 25, 47, 65 | K1-40-2   | 24     |
| A11-50-1    | 69                     | B2-1714-40  | 35, 57, 63         | B4-1212-30   | 65         | K1-40-200 | 24, 46 |
| A11-60-1    | 69                     | B2-188-32   | 67                 | B4-1212-60   | 65         | K1-40-3   | 24     |
| A11-70-1    | 69                     | B2-1913-34  | 67                 | B4-2424-60   | 65         | K1-40-4   | 24     |
| A11-80-1    | 69                     | B2-247-30   | 67                 | B4-33-30     | 65         | K1-40-5   | 24     |
| A11-90-1    | 69                     | B2-2525     | 23, 29, 45, 51, 61 | B4-44-20     | 65         | K1-40-6   | 24     |
| A1-20-1     | 68                     | B2-2525-100 | 23, 45, 61         | B4-66-20     | 65         | K1-40-7   | 24     |
| A1-21-1     | 68                     | B2-2525-30  | 23, 29, 45, 51, 61 | B4-66-40     | 65         | K1-40-8   | 24     |
| A2-10-1     | 68                     | B2-2525-50  | 23, 29, 45, 51, 61 | B4-84-40     | 65         | K1-50-1   | 26     |
| A2-10-2     | 68                     | B2-2525-70  | 23, 45, 61         | B4-88-40     | 65         | K1-50-110 | 26, 48 |
| A2-20-1     | 68                     | B2-2712     | 31,53,63           | B4-88-60     | 65         | K1-50-170 | 26, 48 |
| A2-20-2     | 68                     | B2-2712-100 | 31,53,63           | B5-2710-44-D | 66         | K1-50-2   | 26     |
| A3-10-1     | 69                     | B2-2712-40  | 31,53,63           | B5-278-44-S  | 66         | K1-50-3   | 26     |
| A3-10-2     | 69                     | B2-2717     | 35, 57, 63         | B5-298-55-D  | 66         | K1-50-4   | 26     |
| A5-10-1     | 69                     | B2-2717-100 | 35, 57, 63         | B5-298-55-S  | 66         | K1-50-5   | 26     |
| A5-10-2     | 69                     | B2-2717-40  | 35, 57, 63         | B5-4510-44-D | 66         | K1-50-6   | 26     |
| A5-10-3     | 69                     | B2-2718-32  | 67                 | B5-458-44-S  | 66         | K1-50-7   | 26     |
| A5-10-4     | 69                     | B2-2718-42  | 67                 | B5-468-55-S  | 66         | K1-50-8   | 26     |
| A5-10-5     | 69                     | B2-4025     | 21, 43, 62         | B5-468-55-D  | 66         | K1-60-1   | 32     |
| A5-10-6     | 69                     | B2-4025-100 | 21, 43, 62         | B7-5728-135  | 67         | K1-60-100 | 32, 54 |
| A5-10-7     | 69                     | B2-4025-30  | 21, 43, 62         | B7-5728-265  | 67         | K1-60-130 | 32, 54 |
| A5-10-8     | 69                     | B2-4025-50  | 21, 43, 62         | B8-6030-100  | 67         | K1-60-2   | 32     |
| A6-10-1     | 69                     | B2-4025-70  | 21, 43, 62         | B8-6030-150  | 67         | K1-60-3   | 32     |
| A6-10-2     | 69                     | B2-4825     | 62                 | B8-6030-40   | 67         | K1-60-4   | 32     |
| A7-2525     | 22, 29, 45, 51, 69     | B2-4825-100 | 62                 | B8-6040-100  | 67         | K1-60-40  | 32     |
| A7-2712     | 31, 53, 69             | B2-4825-30  | 62                 | B8-6040-150  | 67         | K1-60-5   | 32     |
| A7-4025     | 21,43,69               | B2-4825-50  | 62                 | B8-6040-200  | 67         | K1-60-6   | 32     |
| A7-5425     | 19, 41, 69             | B2-4825-70  | 62                 | K1-10-1      | 18         | K1-60-7   | 32     |
| B1-1714     | 35, 57, 63             | B2-5333     | 27, 49, 63         | K1-10-110    | 18, 40     | K1-60-70  | 32, 55 |
| B1-1714-100 | 35, 57, 63             | B2-5333-100 | 27, 49, 63         | K1-10-135    | 18, 40     | K1-60-8   | 32     |
| B1-1714-40  | 35, 57, 63             | B2-5333-30  | 27, 49, 63         | K1-10-150    | 18, 40     | K1-70-1   | 32     |
| B1-2525     | 23, 29, 45, 51, 61, 65 | B2-5333-50  | 27, 49, 63         | K1-10-2      | 18         | K1-70-100 | 32, 55 |
| B1-2525-100 | 23, 45, 61             | B2-5333-70  | 27, 49, 63         | K1-10-200    | 18, 40     | K1-70-130 | 32, 55 |
| B1-2525-30  | 23, 29, 45, 51, 61     | B2-5425     | 19,41, 62          | K1-10-260    | 18, 40     | K1-70-2   | 32     |
| B1-2525-50  | 23, 29, 45, 51, 61     | B2-5425-100 | 19,41, 62          | K1-10-3      | 18         | K1-70-3   | 32     |
| B1-2525-70  | 23, 45, 61             | B2-5425-30  | 19,41, 62          | K1-10-4      | 18         | K1-70-4   | 32     |
| B1-2712     | 31, 53, 61             | B2-5425-50  | 19,41, 62          | K1-10-5      | 18         | K1-70-40  | 32     |
| B1-2712-100 | 31, 53, 61             | B2-5425-70  | 19,41, 62          | K1-10-6      | 18         | K1-70-5   | 32,32  |
| B1-2712-40  | 31, 53, 61             | B2-6625     | 25, 47, 63         | K1-10-7      | 18         | K1-70-6   | 32     |
| B1-2717     | 33, 55, 61             | B2-6625-100 | 25, 47, 63         | K1-10-8      | 18         | K1-70-7   | 32     |
| B1-2717-100 | 33, 55, 61             | B2-6625-30  | 25, 47, 63         | K1-20-1      | 20         | K1-70-70  | 32, 55 |
| B1-2717-40  | 33, 55, 61             | B2-6625-50  | 25, 47, 63         | K1-20-110    | 20, 42     | K1-70-8   | 32     |
| B1-4025     | 21, 43, 60, 64         | B2-6625-70  | 25, 47, 63         | K1-20-135    | 20, 42     | K1-80-1   | 34     |
| B1-4025-100 | 21, 43, 60             | B2-84-28    | 67                 | K1-20-150    | 20, 42     | K1-80-100 | 34, 56 |
| B1-4025-30  | 21, 43, 60             | B3-2525     | 23, 45, 65         | K1-20-2      | 20         | K1-80-130 | 34, 56 |
| B1-4025-50  | 21, 43, 60             | B3-2525-100 | 23, 45, 65         | K1-20-200    | 20, 42     | K1-80-2   | 34     |
| B1-4025-70  | 21, 43, 60             | B3-2525-50  | 23, 45, 65         | K1-20-3      | 20         | K1-80-3   | 34     |
| B1-4825     | 60, 64                 | B3-2525-70  | 23, 45, 65         | K1-20-4      | 20         | K1-80-4   | 34     |
| B1-4825-100 | 60                     | B3-4025     | 21, 45, 65         | K1-20-5      | 20         | K1-80-40  | 34     |
| B1-4825-30  | 60                     | B3-4025-100 | 21, 43, 65         | K1-20-6      | 20         | K1-80-5   | 34     |
| B1-4825-50  | 60                     | B3-4025-50  | 21, 43, 65         | K1-20-7      | 20         | K1-80-6   | 34     |
| B1-4825-70  | 60                     | B3-4025-70  | 21, 43, 65         | K1-20-8      | 20         | K1-80-7   | 34     |
| B1-5333     | 27, 49, 61             | B3-4825     | 67                 | K1-30-1      | 22         | K1-80-70  | 34, 56 |
| B1-5333-100 | 27, 49, 61             | B3-4825-100 | 67                 | K1-30-110    | 22, 44     | K1-80-8   | 34     |
| B1-5333-50  | 27, 49, 61             | B3-4825-50  | 67                 | K1-30-135    | 22, 44     | K1-90-1   | 28     |
| B1-5333-70  | 27, 49, 61             | B3-4825-70  | 67                 | K1-30-150    | 22, 44     | K1-90-2   | 28     |
| B1-5425     | 19, 41, 60, 64         | B3-5333     | 27, 49, 65         | K1-30-2      | 22         | K1-90-3   | 28     |
| B1-5425-100 | 19,41, 60              | B3-5333-100 | 27, 49, 65         | K1-30-200    | 22, 44     | K1-90-4   | 28     |
| B1-5425-30  | 19,41, 60              | B3-5333-50  | 27, 49, 65         | K1-30-260    | 22, 44     | K1-90-45  | 28     |
| B1-5425-50  | 19,41, 60              | B3-5333-70  | 27, 49, 65         | K1-30-3      | 22         | K1-90-5   | 28     |
| B1-5425-70  | 19,41, 60              | B3-5425     | 19, 41, 65         | K1-30-4      | 22         | K1-90-6   | 28     |
| B1-6625     | 25, 47, 61, 65         | B3-5425-100 | 19, 41, 65         | K1-30-5      | 22         | K1-90-7   | 28     |

| Code     | Page   | Code       | Page | Code     | Page |
|----------|--------|------------|------|----------|------|
| K1-90-75 | 28, 50 | K3-40-2    | 46   | KS1-10-3 | 18   |
| K1-90-8  | 28     | K3-40-3    | 46   | KS1-10-4 | 18   |
| K2-10-1  | 40     | K3-40-4    | 46   | KS1-10-5 | 18   |
| K2-10-2  | 40     | K3-40-5    | 46   | KS1-10-6 | 18   |
| K2-10-3  | 40     | K3-40-6    | 46   | KS1-10-7 | 18   |
| K2-10-4  | 40     | K3-40-7    | 46   | KS1-10-8 | 18   |
| K2-10-5  | 40     | K3-40-8    | 46   | KS1-20-1 | 20   |
| K2-10-6  | 40     | K3-50-1    | 48   | KS1-20-2 | 20   |
| K2-10-7  | 40     | K3-50-2    | 48   | KS1-20-3 | 20   |
| K2-10-8  | 40     | K3-50-3    | 48   | KS1-20-4 | 20   |
| K2-20-1  | 42     | K3-50-4    | 48   | KS1-20-5 | 20   |
| K2-20-2  | 42     | K3-50-5    | 48   | KS1-20-6 | 20   |
| K2-20-3  | 42     | K3-50-6    | 48   | KS1-20-7 | 20   |
| K2-20-4  | 42     | K3-50-7    | 48   | KS1-20-8 | 20   |
| K2-20-5  | 42     | K3-50-8    | 48   | KS1-30-1 | 22   |
| K2-20-6  | 42     | K3-60-1    | 54   | KS1-30-2 | 22   |
| K2-20-7  | 42     | K3-60-2    | 54   | KS1-30-3 | 22   |
| K2-20-8  | 42     | K3-60-3    | 54   | KS1-30-4 | 22   |
| K2-30-1  | 44     | K3-60-4    | 54   | KS1-30-5 | 22   |
| K2-30-2  | 44     | K3-60-5    | 54   | KS1-30-6 | 22   |
| K2-30-3  | 44     | K3-60-6    | 54   | KS1-30-7 | 22   |
| K2-30-4  | 44     | K3-60-7    | 54   | KS1-30-8 | 22   |
| K2-30-5  | 44     | K3-60-8    | 54   | T1-10-1  | 70   |
| K2-30-6  | 44     | K3-70-1    | 52   | T1-10-2  | 70   |
| K2-30-7  | 44     | K3-70-2    | 52   | T1-10-3  | 70   |
| K2-30-8  | 44     | K3-70-3    | 52   | T1-10-4  | 70   |
| K2-40-1  | 46     | K3-70-4    | 52   | T1-10-5  | 71   |
| K2-40-2  | 46     | K3-70-5    | 52   | T1-10-6  | 71   |
| K2-40-3  | 46     | K3-70-6    | 52   | T1-10-7  | 71   |
| K2-40-4  | 46     | K3-70-7    | 52   | T1-10-8  | 71   |
| K2-40-5  | 46     | K3-70-8    | 52   |          |      |
| K2-40-6  | 46     | K3-80-1    | 56   |          |      |
| K2-40-7  | 46     | K3-80-2    | 56   |          |      |
| K2-40-8  | 46     | K3-80-3    | 56   |          |      |
| K2-50-1  | 48     | K3-80-4    | 56   |          |      |
| K2-50-2  | 48     | K3-80-5    | 56   |          |      |
| K2-50-3  | 48     | K3-80-6    | 56   |          |      |
| K2-50-4  | 48     | K3-80-7    | 56   |          |      |
| K2-50-5  | 48     | K3-80-8    | 56   |          |      |
| K2-50-6  | 48     | K3-90-1    | 50   |          |      |
| K2-50-7  | 48     | K3-90-2    | 50   |          |      |
| K2-50-8  | 48     | K3-90-3    | 50   |          |      |
| K3-10-1  | 40     | K3-90-4    | 50   |          |      |
| K3-10-2  | 40     | K3-90-5    | 50   |          |      |
| K3-10-3  | 40     | K3-90-6    | 50   |          |      |
| K3-10-4  | 40     | K3-90-7    | 50   |          |      |
| K3-10-5  | 40     | K3-90-8    | 50   |          |      |
| K3-10-6  | 40     | KP1-10-110 | 18   |          |      |
| K3-10-7  | 40     | KP1-10-135 | 18   |          |      |
| K3-10-8  | 40     | KP1-10-150 | 18   |          |      |
| K3-20-1  | 42     | KP1-10-200 | 18   |          |      |
| K3-20-2  | 42     | KP1-10-260 | 18   |          |      |
| K3-20-3  | 42     | KP1-20-110 | 20   |          |      |
| K3-20-4  | 42     | KP1-20-135 | 20   |          |      |
| K3-20-5  | 42     | KP1-20-150 | 20   |          |      |
| K3-20-6  | 42     | KP1-20-200 | 20   |          |      |
| K3-20-7  | 42     | KP1-30-110 | 22   |          |      |
| K3-20-8  | 42     | KP1-30-135 | 22   |          |      |
| K3-30-1  | 44     | KP1-30-150 | 22   |          |      |
| K3-30-2  | 44     | KP1-30-200 | 22   |          |      |
| K3-30-3  | 44     | KP1-30-260 | 22   |          |      |
| K3-30-4  | 44     | KP1-40-150 | 24   |          |      |
| K3-30-5  | 44     | KP1-40-200 | 24   |          |      |
| K3-30-6  | 44     | KP1-50-110 | 26   |          |      |
| K3-30-7  | 44     | KP1-50-170 | 26   |          |      |
| K3-30-8  | 44     | KS1-10-1   | 18   |          |      |
| K3-40-1  | 46     | KS1-10-2   | 18   |          |      |





İstiklal Mah. 58. Sokak  
No: 61 PK: 55310  
Tekkeköy/Samsun/Türkiye

+90 362 363 00 01  
victorishealth.com  
info@victorishealth.com

Secure Today  
Elevate Tomorrow



All rights reserved. No part of this catalog may be reproduced, copied,  
or published in any form without the written permission of